VL01-0000000830

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A GLOBAL BIOTECHNOLOGY MARKETING COMPANY

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Nucleic Acid BCT and RNA Complete BCT Letter of Equivalence

Nucleic Acid BCT vs RNA Complete BCT report

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Dear Valued Customer,

This is to notify you that Nucleic Acid BCT™ (10 mL) is comprised of the same glass material, stopper material and reagent as RNA Complete BCT®. The volume of reagent in Nucleic Acid BCT (10 mL) is equivalent to RNA Complete BCT. The performance characteristics of each tube are equivalent.

Nucleic Acid BCT (10 mL) is For Research Use Only. Not for use in diagnostic procedures. Nucleic Acid BCT (10 mL) should only be used for research or the development of new assays.

Kind Regards,

A handwritten signature in black ink, appearing to read "Bill Rock".

Bill Rock
Product Manager

Nucleic Acid BCT™ (10 mL) and RNA Complete BCT® Substantial Equivalence Report

A 10 mL configuration of Nucleic Acid BCT has been developed to provide researchers with approximately twice the amount of plasma as the 5 mL configuration.

This report demonstrates that Nucleic Acid BCT (10 mL) has equivalent performance to RNA Complete BCT as it relates to limiting hemolysis, minimizing changes in plasma volume and maintaining plasma cell-free RNA (cfRNA) and extracellular vesicle (EV) concentration during whole blood sample handling and storage.

Nucleic Acid BCT (10 mL) and RNA Complete BCT are For Research Use Only. Not for use in diagnostic procedures and should only be used for research or in the development of new assays.

The following study and associated data analysis were performed by Streck R&D personnel:

Study:

Venous whole blood from 6 self-declared healthy donors was drawn into Nucleic Acid BCT (10 mL) or RNA Complete BCT and plasma was isolated immediately after draw (Draw time) or after 7 days of ambient temperature storage (Day 7). Plasma samples were screened for hemolysis using the NanoDrop OneC spectrophotometer with the included Hemolysis program. Total recovered plasma weight was measured using a zeroed scale and calculated using the published plasma density of 1.025 g/mL. Nucleic acids were purified from 3 mL of plasma using the QIAamp Circulating Nucleic Acid Kit (QIAGEN®), and concentration was measured with the Qubit miRNA Assay (cfRNA). EVs were assayed using the NanoSight NS300. Plasma was diluted into PBS to a concentration within the specified working range of the instrument. Final concentration was then calculated based on dilution used.

Results:

When hemolysis and plasma volume were monitored at draw and after seven days of ambient temperature whole blood storage, we observed that Nucleic Acid BCT (10 mL) and RNA Complete BCT performed almost identically (Figure 1). When cfRNA were analyzed (Figure 2), we observed that both tubes 1) maintained cfRNA concentrations to draw-time levels after seven days of ambient temperature storage and 2) displayed comparable stabilizing performance. A similar trend was observed for EV concentrations, where Nucleic Acid BCT (10 mL) and RNA Complete BCT performed equivalently (Figure 3).


Figure 1. Nucleic Acid BCT (10 mL) and RNA Complete BCT similarly prevent hemolysis and loss of plasma volume during ambient temperature storage. Qualitative (A) and quantitative (B) analysis of plasma volume loss and hemolysis following seven days of ambient storage in Nucleic Acid BCT (NAC) or RNA Complete BCT (BCT). n=6; ns, not significant by paired T-test. Representative image from a single donor is shown. No significant visual differences were observed between NAC and RNAC BCTs for any of the six donors.

Plasma cfRNA Concentration



Figure 2. Plasma cell-free RNA is maintained to near draw-time levels in blood collected into Nucleic Acid BCT (10 mL) and RNA Complete BCT. Concentration of cell-free RNA in plasma isolated from Nucleic Acid BCT (NAC) or RNA Complete BCT (RNAC) at draw and after seven days of ambient temperature storage. n=6; ns, not significant by paired T-test.

Figure 3. Plasma EV concentrations are maintained to near draw-time levels in blood collected into Nucleic Acid BCT (10 mL) and RNA Complete BCT. Concentration of EVs in plasma isolated from Nucleic Acid BCT (NAC) or RNA Complete BCT (RNAC) at draw and after seven days of ambient temperature storage. n=6; ns, not significant by paired T-test.

Plasma EV Concentration

Conclusion:

Overall, the results of this analysis suggest that Nucleic Acid BCT (10 mL) and RNA Complete BCT have analytically equivalent stabilization performance as it relates to limitation of hemolysis and loss of plasma volume and preservation of draw-time cfRNA and EV concentration.